



Effects far from equilibrium in electromagnetic heating of tissues



Dr. Oliver Szasz^{1,2}, Dr. Gabor Andocs³, Dr. Andras Szasz^{1,2}


(1) Oncotherm GmbH, Troisdorf, Germany

(2) Department of Biotechnics, Faculty of Engineering, St. Istvan University, Budapest, Hungary,

(3) Department of Pharmacology and Toxicology, Faculty of Veterinary Science, St Istvan University, Budapest, Hungary,

Objective

One of the very first treatment "technologies" for oncology is the regional heating of the tissues and body parts, (hyperthermia, HT) This long history was not enough to be accepted as conventional treatment, facing mostly skeptic opinions among the oncology experts. The main reason is its controversial results and poor control, the missing of appropriate selective, controllable safe deep heat delivery. Constrained balance of physiological feedback and the sophisticated transport network with very heterogeneous tissue structures block applying the simple heating practices. This situation requests definite bioengineering tasks as well as new paradigm for the medical applications. The modern heating technologies based on electromagnetic interactions made a huge step ahead in this complex field, but not enough yet to solve some crucial problems in deep heating [1]. The commonly applied microwave and high radio frequency (RF) radiation is challenged by the magnetic and capacitive heating techniques, applied lower frequencies. oncothermia method (OTM), heats by the dielectric loss in the various body electrolytes and makes possible to select between the tissues and concentrate on the malignant cells. Due to the constrained RF current conduction of oncothermia, the complex impedance determines the actual flow direction of the current. This could distinguish by the microscopic heterogeneity in the treated tissue [2]. The main problem is the temperature, which would like to equalized by time in the heated area, and steadily heats up the full environment in wider and wider range, supplying the tumor for growth. We need energy input which can be focused and has no physiologic control. This is which oncothermia had introduced.




By time the temperature equalizes in the area, the heat is conducted away

Oncothermia solves the problem: selectively forces various pathways of apoptosis by electric field

Method

The relatively low frequency RF current dominantly flows in the extracellular electrolyte. The energy absorption creates a temperature gradient through the cellular membrane, which drives non equilibrium processes by constrained heat flow through the membrane, [3]. The ion and mass flows could be well approached by Onsager's theory in the frame of non equilibrium thermodynamic description.



Results

Synthesis of HSP-s additional extracellular and membrane HSP70 appears through the more permeable membrane $T_s + \Delta T$, $\Delta T \approx 0.01 \text{ } ^\circ\text{C}/10 \text{ nm} = 10^{-8} \text{ } ^\circ\text{C}/\text{cm}$, **Temperature-gradient driven processes**
Thermo-electrical current = 150 [pA/ μm^2] (Na⁺ influx), normal = 12 [pA/ μm^2] Na⁺ efflux), drastically decreases the membrane potential, destabilizes the membrane
Thermo-mechanical pressure = 1320 kPa, (electro osmotic effect, rigid tumorcell membrane), water pressure

Heat flow = 1.5 [pW/ μm^2] (at 1 [K/s]), (metabolic heat flow = 0.002 [pW/ μm^2], destroys the ordered membrane
Rectifying effect leads a positive feedback to gain the temperature and the pressure in the membrane
Specific absorption rate of water is high in the membrane (Beta dispersion, ~10 MHz)
Membrane associated apoptotic pathways are activated (E cadherin, beta catenin, p53 expression)

Oncothermia paradigm avoids from high temperature, because:
Temperature heats up the vicinity of the tumor, it can not heat locally focused
Temperature increases the danger of burn of healthy parts in depth (misfocusing, conduction, etc)
Temperature requests the increase of the safety coil up on the skin
The increased surface cooling blocks the temperature sensing in the skin,
The increased surface cooling makes the skin even more isolating, and so the electric burn is more likely
Temperature increases the blood flow in the region, in consequence increases the dissemination
In complementary application with radiotherapy the forced high temperature suppresses the efficacy or blocks at all the effect of radiotherapy
In complementary application with chemotherapy the forced high temperature suppresses the efficacy or blocks at all the chemotherapeutic into the tumor (misconcentration or blood vessel blockage in the tumor)
In complementary application with chemotherapy the forced high temperature increases the cytotoxic side effects in the heads of healthy tissues around by increased chemoreaction rates (vasodilatation in the healthy tissues)
The toxins from the necrotic cells are rapidly transported into the whole body, challenging the anyway low immune status of the patient

In oncothermia the temperature is not a correct dose control, because:
Microstate temperature avoids the natural contra regulation effects
Temperature does not exceed the systemic physiological limit (42 °C)
Tumor "sinks" in is isolated by non temperature dependent way (electric concept)
Focus is to be fixed to the tumor, moves together with the natural body movements (impedance control)
Selection is solved on cellular level suppress the dissemination of the malignant cells
Cellular connections (adherent connections, gap junctions) of malignant cells are established to avoid the further dissemination
Cellular communications (social signal) is established to promote the natural (programmed) cell death for malignant cells
Possibility of the cellular molecular exchange (gap junctions) is reestablished to promote the normal function of the cells
The "master switch" (p53 gene) is activated promoting the natural way of various cell killing pathways
Cell membrane permeability is increased to express the HSP on the outer membrane signaling the cell malignancy for the systemic immune actions
Cell membrane is excited to ignite various communication pathways in the cells
Electric field blocks the positive feedback loop of tumor supporting injury currents


Oncothermia avoids a static equilibrium
Measurement of intensive thermodynamic parameters (like temperature) supposes at least local equilibrium, which never could be reached due to the intensive contra regulatory effects (This concept however, became the main request of the classical hypothermia approach in its guidelines)
The forced equilibrium increases the heat flow to the blood stream, which is an effective cooling media trying to block the static concept
The heat flow to blood supports the positive feedback loop of the basic acidic electrolyte balance, and promotes the intensive growth of the tumor by addressed oxygen delivery
Static constrains try to block the natural dynamics of the living system, which mobilizes its forces to keep the dynamic equilibrium instead of the static one. This creates protection mechanisms of the actual status quo in the tissues, defending the tumor instead its administration (These processes like intracellular HSP development, like forced delivery of metabolic species [oxygen and nutrient], like systemic cooling control, like various stress reactions, etc.)
Process matching equilibrium mobilizes higher level of physiological contra actions and accelerates a competitor between the constants and the nature. This falsely mobilizes the natural healing forces. (Natural actions are gained against the actual treatment and not against the "common enemy", against the malignancy)

Oncothermia works with entirely dynamic (natural) processes
Oncothermia uses tumor killing approach, which is well fitted to the dynamics of the living system, does not constrain for false defence
Control of oncothermia is natural, always fitted to the actual conditions (changes of the electrolytes determines its actions)
No considerable heat flow to the blood stream by oncothermia, no gain of the positive feedback of electrolyte balancing loop
Thermal gradients make dynamism in a very local area of the cell membrane of malignant cells. The applied selection focuses on this thermal non equilibrium
The relatively slow "step up" heating keeps the non equilibrium stable for long time for action
The slow heating up does not create considerable physiological contra actions,
The slow heating makes the healthy tissue adapted to the growing temperature
The slow temperature change does not generate high stress and following stress reactions
The applied electric field makes at least three times more effective cell killing than the temperature does
The applied fractal modulation makes possible selecting and supporting the natural processes to activate the natural healing mechanisms and reestablish the healthy "social signal" between the isolated cells, promoting the anti malignancy on cell level
Complete relaxation could be supported by relaxing music, video or sound effects during the treatment

Oncothermia is far from ideal
Oncothermia is mainly regulated by the patient's tolerance
Oncothermia control based on thermal sensing of the patients, for safety and for efficacy reasons. Safety is avoid burning the tissue of the subcutaneous layers, the efficacy to apply such energy, which does not overload the patient natural defend/protective system
Oncothermia uses natural processes to cure, understanding and using these needs hinking doctors and their understandings
Oncothermia acts of natural physiology regulation, which needs understanding of the processes
Oncothermia needs permanent dynamic approach, follow up as i what is happening during the treatment
Step up heating is the basic treatment approach, which requests permanent care on the process
The effect of the activated natural processes are not acting immediately. To have a control treatment by treatment is essential
The patient's well being during and after the treatment is necessary side of the well conducted protocol
Complete relaxation could be supported by relaxing music, video or sound effects during the treatment

Conclusion

The non equilibrium thermodynamics makes OTM feasible to go over the difficulties of the problem of the selective deep heating. With this new paradigm OTM could be a candidate in the branch of modern therapies in medical practice. OTM could be applied in various biomedical fields where the selection and the drug targeting as well as the personalized treatment are important requests.



References

- [1] Szasz A. Physical background and technical realization of hyperthermia, in Locoregional Radiofrequency Perfusional and Wholebody Hyperthermia in Cancer Treatment. New clinical aspects, (Eds Baronzio GF, Hager ED), Springer Science Eureka.com, Ch 3 pp 27-59, 2006
- [2] Szasz A. Hyperthermia a Modality in the Wings, Journal Cancer Research and Therapeutics, 3:56-66, 2007
- [3] Szasz A, Vincze Gy, Szasz O, Szasz N. An energy analysis of extracorporeal hyperthermia, Magneto and electro biology, 22 (2003) 103-115
- [4] Szasz A, Vincze Gy. Dose concept of oncological hyperthermia heat equation considering the cell destruction, Journal of Cancer Research and Therapeutics, 2:171-181, 2007
- [5] Szasz A, Szasz O, Szasz N. Electrohyperthermia a new paradigm in cancer therapy, Wissenschaft & Forschung, Deutsche Zeitschrift für Onkologie, 2001; 33:91-99
- [6] Szasz A, P. Vincze G, Szasz A. Pink noise behaviour of biosystems; Eur Biophysics J. 30:227-231, 2001
- [7] Vincze Gy, Szasz A, Szasz N. On the thermal noise limit of cellular membranes, Bioelectromagnetics, 28:28-35, 2005,