

## **Modulated electro-hyperthermia, (mEHT) From LAB to clinic**

**Oliver Szasz, PhD<sup>1\*2</sup>, Andras Szasz, PhD<sup>2</sup>**

(1) Oncotherm GmbH, Troisdorf, Germany, [SzaszOliver@oncotherm.de](mailto:SzaszOliver@oncotherm.de)  
(2) Department of Biotechnics, Faculty of Engineering, St. István University, Budaörs,  
Hungary \*[biotech@gek.szie.hu](mailto:biotech@gek.szie.hu)

**Presented at ESHO, Berlin, 2018**

### **Cite this article as:**

Szasz O and Szasz A. (2018): Modulated electro-hyperthermia, (mEHT)  
From LAB to clinic; Oncothermia Journal 23:24-61  
[www.oncothermia-  
journal.com/journal/2018/Modulated\\_electro\\_hyperthermia\\_mEHT.pdf](http://www.oncothermia-journal.com/journal/2018/Modulated_electro_hyperthermia_mEHT.pdf)

# **Modulated electro-hyperthermia, (mEHT) FROM LAB TO CLINIC**


**Oliver Szasz, PhD<sup>1\*,2</sup>, Andras Szasz, PhD<sup>2</sup>**

(1) Oncotherm GmbH, Troisdorf, Germany, [SzaszOliver@oncotherm.de](mailto:SzaszOliver@oncotherm.de)  
(2) Department of Biotechnics, Faculty of Engineering, St. István University, Budaörs, Hungary  
[\\*biotech@gek.szie.hu](mailto:*biotech@gek.szie.hu)

## **Outline**


- The challenge**
- The methods**
- Experimental results**
- Clinical results**
- Conclusions**

## Challenge in oncology – malignancy is systemic



Is the task to simply heat up the tumor?


### Whole-body hyperthermia (WBH, 40-42°C)




Thrall DE, et al. Radiation plus local hyperthermia versus radiation plus the combination of local and whole-body hyperthermia in canine sarcomas. Int J Radiat Oncol Biol Phys. 34(5):1087-96., 1996

## Challenges with whole-body hyperthermia (40-42°C)

### metastatic colorectal cancer      malignant pleural mesothelioma



Hildebrandt B, et al. Whole-body hyperthermia in the scope of von Ardenne's systemic cancer multistep therapy (sCMT) combined with chemotherapy in patients with metastatic colorectal cancer: a phase I/II study; Int J Hyperthermia, 20:317-333, 2004




**Decreases the effect of chemotherapy**

Bakhshandeh A, Wiedemann G, Zabel P, et al. Randomized trial with ICE (ifosfamide, carboplatin, etoposide) plus whole body hyperthermia versus ICE chemotherapy for malignant pleural mesothelioma. Journal of Clinical Oncology, 2004 ASCO Annual Meeting Proceedings (Post-Meeting Edition). 22:14S, 7288. (Supplement), 2004:

### Clinical challenge: Local control ↔ survival time

#### Breast studies

C. C. Vernon, et al., Int. J. Rad. Onc. Biol. Phys., 35:731-744, 1996




#### good local control




#### no significant survival

Sherar, M., et al. Int. J. Rad. Onc. Biol. Phys., 39, 371-380; 1997



## Non-small-cell lung cancer (NSCLC) RR ↔ OS



Kay CS, Choi IB, Jang JY, Choi BO, Kim A, Shinn KS. (1996) Thermoradiotherapy in the Treatment of Locally advanced Non-Small Cell Lung Cancer, *J.Korean Soc. Ther. Radiol. Oncol.* 14:115-122

Initial site of disease progression after treatment			P-value
	RT (n = 40)	RT + HT (n = 40)	
No recurrence	31	4	
Primary tumor and/or regional lymph nodes	15	7	
Distant metastasis	2	10	0.07
Both locoregional and distant*	3	4	
Unknown/missing	17	15	


\* Patients in whom the interval between locoregional disease progression and distant metastasis was less than or equal to 1 month

**Distant metastases are induced by conventional hyperthermia**

## Clinical challenge: Local control ↔ survival time


### Cervix uteri (locally advanced cancer)

Zolciak-Siwinska A, et al; (2013) HDR brachytherapy combined with interstitial hyperthermia in locally advanced cervical cancer patients initially treated with concomitant radiochemotherapy – a phase III study; *Radiotherapy and Oncology* 109:194–199



### Superficial tumors study (RT+HT)

Jones EL, et al. Randomized trial of hyperthermia and radiation for superficial tumors; *J Clin Oncol.* 2005; 23:3079-85



## Challenge of dose of hyperthermia

### **CEM43°CT<sub>x</sub> Calibrated in vitro**

Sapareto and Dewey.CEM43, dose model (Sapareto SA, Dewey WC; (1984) THERMAL DOSE DETERMINATION IN CANCER THERAPYInt.J.Rad. Oncol. Biol. Phys. 10:787-800)

$$CEM43^{\circ}CT_x = \sum_{\text{sessions}} \int_0^{t.\text{time}} R^{(43^{\circ}\text{C}-T_x)} dt$$

Unit: °C  
Unit: ??

### Fits to the clinical data

Francena M, et al: Hyperthermia dose-effect relationship in 420 patients with cervical cancer treated with combined radiotherapy and hyperthermia. Eur. J. Cancer, 45:1969-1978 (2009)

$$TRISE = \frac{1}{\text{sessions}} \sum_{\text{sessions}} \int_0^{t.\text{time}} \frac{(T_{50} - 37^{\circ}\text{C})}{\text{treat. time}} dt$$

duration

Rough average  
(unprecise, incorrect)

corresponding principle

$SAR(t) = c \left( \frac{dT}{dt} \right)$

Actual absorption  
(correct approximation)

### The correct dose

$$\text{Absorbed energy} = \sum_{\text{sessions}} \int_0^{\text{duration}} (SAR(t)) dt$$


measured in Gy (J/kg)  
(like in ionizing radiation)

**Technical requirement:  
high efficacy of energy absorption**

### Serious technical challenge


#### Temperature measurement in a phantom

Fenn AJ, King GA. Adaptive radiofrequency hyperthermia-phased array system for improved cancer therapy: phantom target measurements. Int. J. Hyperthermia. 10:189-208 (1994)



#### Efficacy of heating<10%


Here the measurement of the temperature is the only way to approximate the absorbed energy



#### Consequences

1. The energy can not be used like it is in the case of ionizing radiation
2. The **efficacy** of the energy-absorption must be drastically increased to use the absorbed energy as dose

## Scientific definition of thermal effect



## Outline

### □ The challenge

### □ The methods

---


### □ Experimental results

### □ Clinical results

### □ Conclusions


## Heterogeneous vs. homogeneous

**Conventional homogenic heating**  
(radiation + convection + conduction)



The tumor is rather heterogeneous, only small sizes with invasive heating could be heated approximately homogeneously


**Selective, heterogenous heating**  
(targeted molecules of selected cells)  
**mEHT**



The tumor-cells have definite differences in their bioelectromagnetic/biophysical properties allowing their selected heating

**High efficacy, possibility to use the Gy as dose**

## Abscopal: apoptosis and immunogenic cell-death




Vila M, Przedborski S; Neurological diseases: Targeting programmed cell death in neurodegenerative diseases, *Nature Reviews Neuroscience* 4:365–375 (2003)

Derer A, et al. Radio-immunotherapy-induced immunogenic cancer cells as basis for induction of systemic anti-tumor immune responses – pre-clinical evidence and ongoing clinical applications. *Front. Immunol.* 6:505, (2015)


## Bioelectrodynamical selection

Oncological hyperthermia is a method to destroy malignant cells by heat-inducing absorbed energy



## Bioelectrodynamical energy absorption

Oncological hyperthermia is a method to destroy malignant cells by heat-inducing absorbed energy



## Outline


### □ The challenge

### □ The methods

### □ Experimental results


### □ Clinical results

### □ Conclusions




## Morphological change after the treatment **4 hours**

Andocs G, Kondo T. Toyama University, Japan, (2018)




## Morphological change after the treatment **24 hours**

Andocs G, Kondo T. Toyama University, Japan, (2018)





## mEHT effect of modulation

**Animal model:** Colon26 (murine colorectal cancer) cell line derived orthotopic allograft in Balb/C mouse (double tumor model)  
Andocs G, Kondo T. Toyama University, Japan, (2018)




One of the apoptotic pathways (caspase-dependent)






Tsang Y-W, et al. Improved Immunological Tumor Microenvironment by Combined Electro-Hyperthermia followed by Dendritic Cell Immunotherapy; Cancer Immunology, BMC Cancer 15:708, <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC2647246/>

## Apoptosis



Yang K-L, Huang C-C, Chi M-S, Chiang H-C, Wang Y-S, Andocs G, et.al. (2016) In vitro comparison of conventional hyperthermia and modulated electro-hyperthermia, Oncotarget, doi: 10.18632/oncotarget.11444


## DAMPs: calreticulin



Vancsik T, Kovago Cs, Kiss E, Papp E, Forika G, Benyo Z, Meggyeshazi N, Krenacs T. (2018) Modulated electro-hyperthermia induced loco-regional and systemic tumor destruction in colorectal cancer allografts, J Cancer, 9(1): 41-53

## DAMPs: Comparison of the induced HSP70

### Intracellular level of HSP70 expression




All kinds of hyperthermia could induce intracellular HSP70 expression.


Yang K-L, et al. (2016). Oncotarget, doi: 10.18632/oncotarget.11444

### Extracellular level of HSP70 expression

Yang K-L, Huang C-C, Chi M-S, Chiang H-C, Wang Y-S, Andocs G, et al. (2016) In vitro comparison of conventional hyperthermia and modulated electro-hyperthermia, Oncotarget, doi: 10.18632/oncotarget.11444



### Membrane localization of HSP70 (14h)




Andocs G, et al. (2014). Cell Stress and Chaperones 20(1):37-46

Oncothermia triggered a significantly secretion of HSP70 from cancer cells.


## DAMPs: HMGB1

48 h  
48 h




Vancsik T, Kovago Cs, Kiss E, Papp E, Forika G, Benyo Z, Meggyeshazi N, Krenacs T. (2018) Modulated electro-hyperthermia induced loco-regional and systemic tumor destruction in colorectal cancer allografts, J Cancer, 9(1): 41-53

## DAMP-process: Antigen presenting cells




Vancsik T, Kovago Cs, Kiss E, Papp E, Forika G, Benyo Z, Meggyeshazi N, Krenacs T. (2018) Modulated electro-hyperthermia induced loco-regional and systemic tumor destruction in colorectal cancer allografts, *J Cancer*, 9(1): 41-53




Marsdenia tenacissima (MTE) injection before mEHT

Chlorogenic-acid (11 mg/ml)



### Abscopal effect by Marsdenia tenacissima




Vancsik T, et al (2018) *Journal of Cancer* 9:41-53.

## mEHT heats the cell-membrane rafts


**U937**

(human lymphoma cell-line)

**in-vitro**





**HT29, (human CRC cell-line)**



**HT29, (human CRC cell-line); in vitro**

Vancsik T et al., 33rd annual Conference of International Clinical Hyperthermia Society (ICHS), Nidda, Germany, July 11-13, 2015

### temperature sensing dye on cell membrane



## Oncotherapy is thermal


Post-treatment 24 H, 30 min; Each point is 3 independent measurements

**In vitro:** U937 human lymphoma cell-line ( $10^6$  cell/ml)




Andocs G, et al, Biology and Medicine 7(4):1-9, 2017,

**In vivo:** HT29 human colorectal cancer xenograft


Andocs G et al. Radiology and Oncology (Strahlentherapie und Onkologie) 185:120-26 , (2009)



### Arrhenius plots (thermal effect)




Power-function



The activation energy is calculated

## Block the invasion and dissemination



Andocs G, Szasz O, Szasz A (2009) Oncothermia treatment of cancer: from the laboratory to clinic. Electromagn Biol Med 28(2):148–165

Yang K-L, Huang C-C, Chi M-S, Chiang H-C, Wang Y-S, Andocs G, et.al. (2016) In vitro comparison of conventional hyperthermia and modulated electro-hyperthermia, Oncotarget, doi: 10.18632/oncotarget.11444

## Outline



- The challenge
- The methods
- Experimental results
- Clinical results
- - Conclusions

## mEHT treatment – easy to use and safe



### Triple negative locally advanced breast cancer

**Investigator:** Prof. Chi K-W, Shih-Kong Hospital, Taipei, Taiwan  
**Presented on** 35<sup>th</sup> ICHS Conference, Guangzhou, China; Nov.2017)




## Brain metastasis from breast cancer

**Investigator:** Dr. Marwan Akasheh; **Institute:** Dar Alshefa' Tumors Treatment Center, Amman, Jordan, **Patient:** female 53 y.  
**mEHT Monotherapy**



## NSCLC Cases


**Investigator:** Prof.Dr.Chang Geol Lee; **Institute:** Department of Radiation Oncology, Yonsei Cancer Center, Yonsei University Health System, Seoul, South Korea; **Presented on:** KOSG'2013, Jeju  
**Patient:** 51 years old, male patient; **Diagnosis:** NSCLC, RUL, AdenoCa, cT4N3M0; **T:** 6.5x3.7cm, right brachial plexus invasion; **N:** Both SCL, Lt. axillary LN; **Referred for:** CCRT  
**Treatment:** RT: 66 Gy/30 fractions; CTx: weekly paclitaxel/cisplatin; Oncothermia: 10 fractions/5 weeks



## Recurrent uterine sarcoma with peritoneal seedings

Investigator: Prof. Chi K-W, Shih-Kong Hospital, Taipei, Taiwan  
Presented on 35<sup>th</sup> ICHS Conference, Guangzhou, China; Nov.2017)


refractory to chemotherapy and salvage  
with combined radiotherapy (45Gy/30fx)



Intratumoral ipilimumab 2.5 mg, i.v. nivolumab 50 mg and complementary with oncothermia 6 times (1 time/week)

## Recurrent nasopharyngeal carcinoma

Investigator: Prof. Chi K-W, Shih-Kong Hospital, Taipei, Taiwan  
Presented on 35<sup>th</sup> ICHS Conference, Guangzhou, China; Nov.2017




Nivolumab 60mg q3w; RT: 2017/05/01- 6/12:60Gy/30fx; Oncothermia 6x

## Abscopal effect

Investigator: YH Kim; Ewha Womans University Mokdong Hospital, Seoul, Korea



### Recurrent refractive progressive ovarian cancer (55y)




### Metastatic non-small-cell lung cancer (55y)

Investigator: Prof. Dr. Seong Min Yoon,  
Division of Hematology-Oncology, Department of Internal Medicine,  
Samsung Changwon Hospital, Sungkyunkwan University, Korea  
Patient: SAsc, 72 y, male, Primer-tumor: Non-small cell lung cancer;  
Size: 9.5 cm right middle lobe; Metastases: in sentinel and distant  
lymph-nodes; Tumor-classification: cT2 cN2 Mx, stage IIIIB  
Treatment: 28x1.7 Gy; support: **250 microgram Leukine and**  
**Oncothermia 6x**


### Invasive adenocarcinoma of ovary (grade 2) (33 y). Vaginal bleeding; G5P2




Randomized study (n=6+6) for pharmacokinetics



Lee SY, Kim M-G (2015); Int J Hyp, 31:869; 2015



Despite the mild overall heating mEHT significantly increases the maximum of the kinetic curve ( $C_{\text{max}}$ ) and also significantly decreases the time at maximum ( $T_{\text{max}}$ ) parameters. The area under the peak (AUC) was stable, indicating the unchanged systemic adverse effects, despite the increase of the absorption of the drug.




This effect depends on the permeability of the vessel-wall, which is increased by the applied electric field.

## Prospective study (n=20) for blood-flow changes in cervix

Biopsy proven cervical carcinoma by mEHT therapy. **Peri-tumortemperature:** internal organ temperature probe. **Tumor blood flow:** by determining the peak systolic velocity/end-diastolic velocity ratio (S/D ratio) and the resistance index (RI) within blood vessels. (3D color Doppler ultrasound)



Lee S-Y, Kim JH, Han YH, Cho DH; Int. J. Hyperthermia, Jan 3:1-24, (2018)



D. Fatehi, *Technical Quality of Deep Hyperthermia Using the BSD-2000*. Uitgeverij Box Press, Oisterwijk, The Netherlands, 2007, PhD.

## Recurrent glioblastoma multiforme meta-analysis

Roussakow S Clinical and economic evaluation of dose-dense temozolomide 21/28d regimen with and without concurrent modulated electro-hyperthermia in the treatment of recurrent glioblastoma: a retrospective comparison of cohort trials (2a level evidence) , [based on Publications of Gronemeyer et al. (2004)]; BMJ Open, 7:e017387.doi:1136/BMJ-open-2017-017387; (2017).




No Study, cohort	M/No	Mean survival time	Weight (%)							p-value	
			1	2	3	4	5	6	7		
1 Jungk (2016), BCNU	9,0/34	9,00 (7,49 – 10,51)			10,6%	1,00	0,00	0,87	0,98	0,61	0,07
2 Reithmeier (2010), BCNU	5,1/35	5,10 (4,48 – 5,72)			62,6%	0,00	1,00	0,14	0,03	0,00	0,05
3 Glas (2009), ACNU + VM26	6,0/35	6,00 (5,01 – 6,99)			24,4%	0,00	0,14	1,00	0,12	0,02	0,00
4 Heiland (2016), BEV + CCNU	8,7/18	8,70 (5,51 – 11,89)			2,4%	0,87	0,03	0,12	1,00	0,90	0,62
5 Hau (2010), ACNU+mEHT, salvage	9,0/15	8,96 (6,69 – 11,23)			4,7%	0,98	0,00	0,02	0,90	1,00	0,65
6 Douwes (2005), ACNU+mEHT	9,7/19	9,68 (7,61 – 11,75)			5,6%	0,61	0,00	0,62	0,65	1,00	0,04

Random effect model ( $I^2=0,0\%$ ,  $p=0,59$ )



## Glioblastoma multiform – comparison of three survival results

### Overall survivals




### Glioblastoma multiform


**Investigator:** Dr. Gurdev Parmar

**Institute:** Integrated Health Clinic, Cancer care center, Fort Langley, British Columbia, Canada.

**Published:** 33<sup>rd</sup> ICHS Nidda, Germany; 2015



## Glioblastoma multiform – comparison of three survival results



## Small-cell-lung-cancer (n=9+10) double arm prospective study

**Investigator:** Professor DY Lee, Kagnam Severance Hospital, Yonsei University, Seoul, S.Korea


**Published:** Lee DY, et al. (2013) Conference Papers in Medicine, Vol.2013, Article ID 910363, pp.1-7

Prospective, monocenter, cohort double-arm study of chemotherapy with and without complementary oncotherapy

**Chemotherapy 1<sup>st</sup> line (n=28):** Irinotecan (60 mg/m<sup>2</sup>), Cisplatin (60 mg/m<sup>2</sup>) three times.

**Chemotherapy 2<sup>nd</sup> line (n=19):** Etoposide, (110 mg/m<sup>2</sup>) Cisplatin (70 mg/m<sup>2</sup>)

**Additional oncotherapy in 2<sup>nd</sup> line combination (n=9):** 150 Watt, 1,490.5 kJ, 60 min, every second day, with rise in temperature to 38.5°C–42.5°C. Electrode 30 cm diameter at least 12 sessions were in 1 cycle.




**Both the local control and the overall survival are improved**

## Recurrent cervix double arm (n=20+18), randomized study

Lee SY, Lee NR, Cho D-H, Kim JS; Oncology Letters, <https://doi.org/10.3892/ol.2017.6117>, (2017)

Patients received conventional chemotherapy alone (n=20) compared to the combination to mEHT (n=18). Every patient had chemotherapy [paclitaxel + cisplatin (n=14), paclitaxel + carboplatin (n=10), cisplatin + 5-fluorouracil (n=12), cisplatin alone (n=2)]. Radiotherapy was not possible in this cohort.



**Both the local control and the overall survival are improved**

Phase III randomised cervix trial (n=236) of mEHT with RT  
(interim results (n=160), follow-up is ongoing)

**Investigators:** Minnaar CA; Kotzen JA; Baeyens A. Charlotte Maxeke Johannesburg Academic Hospital, S.Afrika. **Aim:** to enrol 236 participants with FIGO stage IIIB (initial distal parametrium involvement) to IIIB cervical cancer

Statistics	n	%
HIV positive	120	51%
Stage III	157	66.6%

Radiation: 25x2Gy external and 3x8Gy brachytherapy



Chemotherapy: 3x 80mg/m<sup>2</sup> Cisplatin

mEHT (oncothermia): 2x 55min/week (4 weeks)


### Local control

Until now: 6 month Local Disease Control  
160 patients completed 6 month PET scan.

Measured	Radio-chemotherapy				Gain by mEHT (%)	
	with mEHT		without mEHT			
	n	%	n	%		
Complete response	33	47%	27	32%	15%	
6 months survival (n=160)	70	91%	90	81%	10%	
24 months survival (n=114)	55	78%	59	65%	13%	
Progression Free survival	53	76%	55	61%	15%	



Interim report by HIV infection (subgroups)




**Until now, both the local control and the overall survival are improved**

## Pancreas Phase II study (n=26)

Metastatic pancreatic cancer with progressive disease after gemcitabine treatment. In the 2<sup>nd</sup> line the patients received gemcitabine 1000mg/m<sup>2</sup> IV and oxaliplatin 100mg/m<sup>2</sup> IV day 1 (GEMOX) combined with mEHT days 1, 3 and 5 all repeated at 14 days


Characteristics	Enrolled (n=17)
Male	9
Female	8
ECOG Performance status	
ECOG 1	5
ECOG2	12
Stage at study entry	
Liver metastasis	6
Lung metastasis	4
Lymph node metastasis	6
Peritoneal carcinosis	4
Bone metastasis	6
Ascites/pleural effusion	8
Nr. of prior chemotherapy cycles (GEM) - median	5.4
Histopathologic types	
Duct cell carcinoma	11
Acinic cell carcinoma	1
Papillary mucinous carcinoma	2
Signet ring carcinoma	1
Adenosquamous carcinoma	1
Undifferentiated carcinoma	1
Prior regional therapy	
Surgery	6
Radiotherapy	3



Volovat et al.; (2014); Romanian Reports in Physics, 66:166–174

## Advanced, metastatic pancreas Phase II study (n=99)


Phase II retrospective, two centers (A & B), single-arm clinical trial (n=99) for advanced pancreas cancer treated by oncotherapy [1]. Most of the patients had distant metastases (77, 77.8%; A:23.88.7%, B:54.74%) and 40% had multiple metastases. The trial includes a cohort of heavily pretreated patients (3+ lines), and due to the refractory or another fail of the conventional therapies, in this study mEHT was applied as monotherapy. The first and subsequent year survivals were: 1st:50.5%, 2nd: 27.3%, 3rd:15.2%, 4th:8.1%, 5th:3%. These values are significantly higher than the values from the large databases (SEER and Eurocare). The center B had the historical arm of conventional therapies (and palliation) of the cohort, with median overall survival 6.5 m, while the median in active study arm was 12.7 m.



Dani A, et al. (2008) Forum Hyperthermie 1:13–20

## Comparison of pancreas studies

### Metastatic pancreas CA 1y survival [%]




### Advanced, high risk recurrent sarcoma Phase II study (n=24)

After recurrence of first line chemotherapy with doxorubicin new line of chemotherapy (ifosfamide 3000mg/m<sup>2</sup>, day 1–3) and mEHT (1 hour application with temperature between 41.5°C and 42°C, 3 days/week).

The response 88% (partial response 44% patients for 4 m; stable disease 44% patients for 4 m and 5% only 1 m).


Characteristic	Nr. of patients
Performance status	
ECOG 2	4
ECOG 3	14
Site of metastasis	
Lung	8
Liver	11
Bone	7
Hystopathologic Type	
Fibrosarcoma	5
Mixofibrosarcoma	2
Synovial sarcoma	3
Leiomysarcoma	3
Epithelioid Sarcoma	2
Angiosarcoma	3

Volovat et al.; (2014) The results of combination of ifosfamid and locoregional hyperthermia (ehy 2000) in patients with advanced abdominal soft-tissue sarcoma after relapse of first line chemotherapy, *Romanian Reports in Physics*, 66:175–181



## Hepatocellular carcinoma Phase II study (n=21)

A mono-institutional uncontrolled phase II trial was conducted on advanced HCC patients. Treatment was continued until disease progression (PD) or unacceptable drug-related toxicities. Sorafenib treatment interruptions and dose reductions (initially 200 mg twice daily, then reduced to 200 mg once daily) were allowed for drug-related toxicity.



## Side effects, contraindications of oncotherapy

Applicable in most of the cases when gold-standards and/or conventional heat-therapies are not obtainable


Oncotherapy does not gain (mostly suppresses) the side-effects of other complementary applied therapies

### Side effects

- Light erythematous redness (<8%)
- Surface burn (very rare)
- Adipose burn (<3%)

**Sorafenib and locoregional deep electro-hyperthermia in advanced hepatocellular carcinoma.** Adverse events (fatigue, diarrhea), are commonly associated with Sorafenib.

Gadaleta-Calderola G. et al Oncol Lett, 2014 Oct,8(4):1783-1787



Pang CLK et al; Clinicaltrial.gov, Identifier: NCT02638051, 2015, Molecular and Clinical Oncology, 6:723-732, 2017

### Contraindications

- Organ transplants or vanished immune-reactions
- Missing thermal sensing
- Inability to communicate
- Woman in pregnancy
- Pacemaker or other electric implants (?)

Sedation, analgesia or anesthesia is not allowed when oncotherapy is applied!

Side effects	mEHT side effects (glioma treatment)	Rel. val.
Short term asthenia after treatment (< 2 h)	9 %	
Local redness (rubor) of the skin	8 %	
Subcutaneous fibrosis of fatty tissue	1 %	
Skin burn (diameter < 1.5 cm) Grade I-II	2 %	
Headache and vomiting (< 2 h)	12 %	

Sahinbas H. et al. Deutsche Zeitschrift fuer Onkologie 39:154–160, 2007


### High care is necessary

- Epileptic or electric sensitivity
- Ascites or other free-electrolytes
- Large metallic implants

## Outline

- The challenge
- The methods
- Experimental results
- Clinical results
- Conclusions

### mEHT results/hot-topics



**Selectively destroys the malignant cells with high metabolic rate**  
Szasz A, et al. Magneto- and electro-biology 22(2):103–115, (2003)  
**Induces massive apoptosis (produces apoptotic bodies)**  
Meggyeshazi N, et al. Strahlenther Onkol 190:815–822, (2014)  
**Apoptotic signal starts from the cell-membrane**  
Andocs G, et al. Biology and Medicine 7(4):1-9, (2015)  
**Expression of: mHSP70, mHSP90, Trail-R2, CRT, Cyt-C, HMGB1**  
Andocs G, et al. Cell Stress and Chaperones 20(1):37-46, (2014)  
**Regulates different genes than the same temperature heating**  
Andocs G, et al. Cell Death Discovery (Nature Publishing Group), 2, 16039, (2016)  
**Produces immunogenic cell-death (ICD)**  
Vancsik T, et al. J. of Cancer, 20(1):37-46, (2018)  
**Forms APC with ICD process**  
Meggyeshazi N, et al. Hindawi, Conference Papers in Medicine, Volume 2013, Article ID 187835, (2013)

**Reestablishes the cellular connections (E-cadherin)**  
Yang K-L, et al. Oncotarget, doi: 10.18632/oncotarget.11444, (2016)  
**Forms bonding connections ( $\beta$ -catenin)**  
Szasz A Thermal Med 29(1):1-23, (2013)  
**The primary tumor is enveloped by lymphocytes**  
Szasz A, et al. A chapter in book Ed: Huijgol N, ISBN 980-953-307-019-8, InTech., (2013)  
**Activates the neutrophils in the envelop**  
Szasz A, et al. Springer Science, Heidelberg, (2010)


**Creates abscopal effect**  
Qin W, et al. Oncol Rep 32(6):2373-2379, (2014)  
**The rechallenge could not be forced**  
Yuk-Wah Tsang, et al. BMC Cancer 15:708, (2015)  
**Immune (and/or small-vehicles) actions**  
Kleef R, et al. Case Rep Oncol 5:212-215, (2012); Schirrmacher V, et al. Oncology Letters 8:2403-2406, (2014); Schirrmacher V, et al. Immunotherapy 7: 855–860; (2015)  
**Tumor-specific immune-reactions (vaccination)**  
Patent (EU) <http://www.google.com/patents/EP2703001A1?cl=en>,  
Patent (USA) <http://www.freepatentsonline.com/20150217099.pdf>

mEHT gains the local control and the survival simultaneously

#### Clinical proofs for primary and metastatic malignancies (selected 21 from 100+ clinical)

<b>Brain gliomas</b>	1. Wismeth C, et al. J Neurooncol 98(3):395–405, (2010) 2. Fiorentini G, et al In Vivo 20(6A):721–724, (2006) 3. Roussakow S, BMJ Open, 7:e017387, (2017) 4. Hager ED, et al. Deutsche Zeitschrift für Onkologie 38(3):133-135, (2006) 5. Hager ED et al. J Clin Oncol, Ann.Meet.Proc.(ASCO Post-Meet.Ed.) 26:2047, (2008) 6. Sahinbas H, et al. Deutsche Zeitschrift fuer Onkologie 39:154–160, (2007)
<b>Uterus cervix cancer</b>	7. Lee S-Y, et al. Oncology Letters, DOI: 10.3892/ol.2017.6117(2017 ) 8. C Minnaar, et al. Physica Medica 32(2):151-152, (2016)
<b>Pancreas cancer</b>	9. Volovat C, et al. Romanian Reports in Physics 66(1):166-174, (2014) 10. Douwes FR, Biologische Medizin 35:126–130, (2006)
<b>Lung cancer</b>	11. Ou J, et al. European J Pharmaceutical Sciences, 109:412-418, (2017) 12. Lee DY, et al. Advances in Lung Cancer 4:1-7, (2015) 13. Lee DY, et al. Conference Papers in Medicine, Volume 2013, Article ID 910363, (2013)
<b>Liver cancer</b>	14. Hager ED. Deutsche Zeitschrift für Onkologie, 36:132-134, (2004) 15. Gadaleta-Caldarola G, et al. Oncol Lett, 2014 Oct;8(4):1783-1787, (2014)
<b>Sarcoma</b>	16. Jeung TS, et al. Case Reports in Clinical Medicine 4:157-168, (2015) 17. Volovat C, et al. Romanian Reports in Physics, 66(1):175-181, (2014)
<b>Intraperitoneum</b>	18. Pang CLK, et al. Molecular and Clinical Oncology, 6:723-732 (2017)
<b>Physiology</b>	19. Lee S-Y, et al. Int J Hyp, 28:1-6'(2015) 20. Lee S-Y, et al. Int. J. Hyperthermia, online: 21 January, (2018)
<b>Monotherapy</b>	21. Jeung TS, et al. Conf. Papers in Medicine, Vol. 2013, Article ID 392480, (2013)

mEHT fits well to the modern oncology



Thank you for your attention!

## Acknowledgement

We would like to say thank you for the continuous help and support of the following institutes and people:

Semmelweis University, Hungary  
Toyama University, Japan  
Department of Radiation Oncology, Taiwan  
Ewha Womans University Hospital, Korea  
University of the Witwatersrand, Johannesburg, South Africa, ...

Dr. Gabor Andocs  
Dr. Tibor Krenacs  
Dr. Marwan Akasheh  
Prof.Dr.Chang Geol Lee  
Dr. Carrie Minnaar, ...

... and thank you for the many additional institutions and people who has helped us.

This project is supported by the National Research, Development and Innovation Office. Project number: NVKP\_16-1-2016-0042